B.Well®

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RUS

KZ

Automatic Blood Pressure Monitor WA-55 Instruction Manual

Table of contents

1. Introduction	4
1.1. Features of the WA-55.....	4
1.2. Important information about self-measurement.....	4
2. Important information on the subject of blood-pressure and its measurement	5
2.1. How does high blood-pressure arise?.....	5
2.2. Table for classifying blood-pressure values	6
2.3. 3Check technology	7
2.4. Pulse Arrhythmia Detection.....	8
2.5. Traffic Light Indication in the Display.....	9
3. The various components of the blood-pressure monitor	10
4. Putting the blood-pressure monitor into operation	12
4.1. Inserting the batteries and using a mains adapter.....	12
4.2. Tube connection.....	13
4.3. Reading the set date.....	13
4.4. Settings the time and date.....	13
4.5. Specifying user.....	14
5. Carrying out a measurement	14
5.1. Before the measurement.....	14
5.2. Common sources of error.....	14
5.3. Fitting the cuff.....	15
5.4. Measurement mode selection.....	15
5.5. Measuring procedure.....	16
5.6. Discontinuing a measurement	17
5.7. Memory – storage and recall of the measurements.....	17
5.8. Clear all values.....	17
6. Error messages/malfunctions	18
7. Care and maintenance, recalibration	20
8. Reference to standards	21
9. Technical specifications	21
10. Guarantee	22

1. Introduction

1.1. Features of the WA-55

Thank you for purchasing the B. Well upper arm blood pressure monitor WA-55. Designed for convenience and easy operation, this device provides fast and reliable measurement of systolic and diastolic blood pressure as well as heart rate using the oscillometric measurement method.

Your blood pressure is an important parameter that can be used to monitor your health. High blood pressure (hypertension) is a serious health problem that is frequently found in the modern world. This device enables you to monitor your blood pressure regularly, and maintain a record of your blood pressure measurements. You can then use this record to help your physician diagnose and prevent any health problems in the future.

Functional features

- The WA-55 is a fully automatic, digital, dual user, upper arm blood pressure measuring device.
- The unique 3Check technology allows You to measure Your blood pressure and obtain an accurate result even during arrhythmia
- Up-to-date Fuzzy logic technology provides the individual level of cuff inflation for quick, precise and painless result.
- The Pulse Arrhythmia Detection technology with an audible sound signal which warns about the malfunctions of normal rate and periodicity of systole during the measurement
- Big 3-line LCD screen with coloured blood pressure indication scale according to European Society of Hypertension (ESH)
- LCD screen and button backlight
- EA- Cuff for arm, washable
- Time/date functions
- Blood pressure level imaging according to European Society of Hypertension (ESH)
- Two separate memory sections by 60 blood pressure readings for 2 users with a date and time
- Battery life indicator
- Cut-off power independent measurement memory
- Automatic switch off
- Mains adapter port

This device is easy to use and has been proven in clinical studies to provide excellent accuracy.

Before using the WA-55, read this instruction manual carefully and keep it in a safe place. For further questions on the subject of blood-pressure and its measurement, please contact your doctor.

1.2. Important information about self-measurement

- Do not forget: self-measurement means control, not diagnosis or treatment. Unusual values must always be discussed with your doctor. Under no circum-

stances should you alter the dosages of any drugs prescribed by your doctor.

- The pulse display is not suitable for checking the frequency of heart pacemakers!
- In cases of cardiac irregularity (Arrhythmia), measurements made with this instrument should only be evaluated after consultation with the doctor.

Electromagnetic interference:

The device contains sensitive electronic components (Microcomputer). Therefore, avoid strong electrical or electromagnetic fields in the direct vicinity of the device (e.g. mobile telephones, microwave cookers). These can lead to temporary impairment of the measuring accuracy.

2. Important information on the subject of blood-pressure and its measurement

2.1. How does high blood-pressure arise?

The level of blood-pressure is determined in a part of the brain, the so-called circulatory centre, and adapted to the respective situation by way of feedback via the nervous system.

To adjust the blood-pressure, the strength and frequency of the heart (Pulse), as well as the width of circulatory blood vessels is altered. The latter is effected by way of fine muscles in the blood vessel walls.

The level of arterial blood-pressure changes periodically during the heart activity: During the «blood ejection» (Systole) the value is maximal (systolic blood-pressure value), at the end of the heart's «rest period» (Diastole) minimal (diastolic blood-pressure value).

The blood-pressure values must lie within certain normal ranges in order to prevent particular diseases.

2.2. Table for classifying blood-pressure values

Please get acquainted with the table for classifying blood-pressure values (unit: mmHg) according to European Society of Hypertension (ESH). For more information see point 2.5

Table for classifying blood-pressure values (unit: mmHg) according to European Society of Hypertension (ESH).

Range	Systolic blood pressure	Diastolic blood pressure	Measures
Grade 3: severe hypertension	Higher or equal to 180	Higher or equal to 110	Urgently seek medical advice!
Grade 2: moderate hypertension	160-179	100-109	Consult your doctor immediately
Grade 1: mild hypertension	140-159	90-99	Consult your doctor
High normal	130-139	85-89	Consult your doctor
Normal	Lower than 130	Lower than 85	Self-check
Optimal	Lower than 120	Lower than 80	Self-check

► Further information

- If your values are mostly standard under resting conditions but exceptionally high under conditions of physical or psychological stress, it is possible that you are suffering from so-called «labile hypertension». Please consult your doctor if you suspect that this might be the case.
 - Correctly measured diastolic blood-pressure values above 110mmHg require immediate medical treatment.
- There are measures which you can take, not only for reducing a medically established high blood pressure, but also for prevention. These measures are part of your general way of life:

A) Eating habits

- Strive for a normal weight corresponding to your age. Reduce overweight!
- Avoid excessive consumption of common salt.
- Avoid fatty foods.

B) Previous illnesses

- Follow consistently any medical instructions for treating previous illness such as:
- Diabetes (Diabetes mellitus)
- Fat metabolism disorder
- Gout

C) Habits

- Give up smoking completely
- Drink only moderate amounts of alcohol
- Restrict your caffeine consumption (Coffee)

D) Physical constitution:

- After a preliminary medical examination, do regular sport.
- Choose sports which require stamina and avoid those which require strength.
- Avoid reaching the limit of your performance.
- With previous illnesses and/or an age of over 40 years, please consult your doctor before beginning your sporting activities. He will advise you regarding the type and extent of types of sport that are possible for you.

2.3.3 Check technology

- 3 Check technology is a new concept for optimum reliability in self-measurement of blood pressure.
- An advanced measurement accuracy is achieved by the automatic analysis of three successive measurements (see 5.5 Measuring Procedure).
- The new system provides reliable values for the doctor and can be used as the basis for reliable diagnostics and medication therapy for high blood pressure.

Why 3Check?

• Human blood pressure is not stable. The unique 3Check technology is more accurate than the technology of calculation of average result which is often used in the devices of other companies. 3Check technology allows You to measure Your blood pressure and obtain an accurate result even during arrhythmia.

A) Key advantages

Reduction of:

- Insufficient rest prior to measurement
- Movement effects
- Agitation

B) Medical benefits

- Improved accuracy
- Reliable patient self-measurement data for the doctor
- Safe hypertension diagnostic
- Reliable therapy control

C) Measurement sequence

- With the full measurement cycles, the total measurement time remains less than 3 minutes, compared to a single measurement time of 1.5 minute.
- Single results are not displayed.
- Due to the Data Analysis result, a 4th measurement may be applied.

2.4. Pulse Arrhythmia Detection

Appearance of the Arrhythmia indicator

The appearance of the symbol  signifies that a certain pulse irregularity was detected during the measurement. The result can vary from your normal blood pressure. As a rule this is not a cause for concern; however, if the symbol  appears more frequently (e.g. several times per week on measurements performed daily) or if it suddenly appears more often than usual, we recommend you inform your doctor. Arrhythmia during the measurement will not affect it if you carry on measurement in 3Check mode (see point 2.3. for information)

Please show your doctor the following explanation:

The device is an oscillometric blood pressure measuring device that also analyses the pulse frequency during measurement. The accuracy of this device has been clinically validated.


If pulse irregularities occur during measurement, the arrhythmia symbol  is displayed after the measurement. If the symbol appears **more frequently** (e.g. several times per week on measurements performed daily) or **if it suddenly appears more often than usual**, we recommend the patient to seek medical advice.

The device does not replace a cardiac examination, but serves to detect pulse irregularities at an early stage.

2.5. Traffic Light Indication in the Display

The coloured bars on the left-hand edge of the display show you the range within which the indicated blood pressure values lies. Depending on the height of the bar, the readout value is either within the normal (green), borderline (yellow and orange) or danger (red) range. The classification corresponds to the 6 ranges in the Table as defined by the ESH and described on the table of the point 2.2. The recommendations of the European Society of Hypertension (ESH) allow to diagnose and treat the hypertension more effectively and do not contradict World Health Organization recommendations.



Additionally, the colours on the display differ according to the readings. If, for example, the readings are in the range of 1, 2, 3 the display light is green, in the range of 4, 5 the light is yellow, and in the range of 6 the display light is red.


3. The various elements of your blood pressure monitor

The illustration shows the blood-pressure monitor WA-55, consisting of:

a) Measuring unit

Item	Function
1 LCD screen	Your recorded blood pressure, heart rate and time are displayed here.
2 Start/Stop button 	Press to start or stop recording your blood pressure.
3 Memory button 	Press to view the previously recorded data.
4 Time button 	Press to set or view current date and time or to change the user.
5 Arm cuff connector 	Connect the upper arm cuff here and place the cuff around your arm to measure your blood pressure.
6 DC 6V 	Connect to adapter.
7 3Check slide switch 	Select 3Check or single mode measurement here.
8 Battery compartment	Insert the four AA batteries here.


b) Cuff: optional

Cuff: optional


EA- Cuff M – L, for arm circumference 22-42 cm, washable

or

EA-Cuff M, for arm circumference 22-32 cm, washable

c) LCD screen scheme with symbol designation


4. Putting the blood-pressure monitor into operation

4.1. Inserting the batteries and using a mains adapter

Follow these steps to insert four AA batteries in the device.

1. Open the battery compartment cover in the direction shown.
2. Insert the four AA batteries with the correct polarity as indicated.
3. Replace the battery compartment cover.

⚠ Attention!

- If the battery  warning appears in the display, the batteries are flat and must be replaced by new ones.
- After the battery  warning appears, the device is blocked until the batteries have been replaced.
- Please use «AA» Long-Life or Alkaline 1.5V Batteries.
- If the blood-pressure monitor is left unused for long periods, please remove the batteries from the device.


Functional check:

Press and hold the   button to test all the display elements. When functioning correctly all segments appear.

Using a mains adapter

- 1) Plug the mains adaptor into a 100-240V, 50/60Hz power socket.
- 2) Plug the DIN plug into the socket at the right side of the instrument.

No power is taken from the batteries while the mains adaptor is connected to the instrument.


⚠ Warning!

If you need mains adapter, You may purchase it separately. Use only the mains adaptor AD-155. The use of any other adapters can make your warranty void.

4.2. Tube connection

Insert the cuff tube into the opening provided on the left side of the instrument, as shown in the diagram.

**4.3. Reading the set date**

Please press the **⌚** («TIME») button.

4.4. Setting the time and date


This blood-pressure monitor incorporates an integrated clock with date display. This has the advantage, that at each measurement procedure, not only the blood-pressure values are stored, but also the exact moment of the measurement.

After new batteries have been inserted, you must enter the date and current time. For this, please proceed as follows:

- 1) Press the **⌚** («TIME») button and hold for 3 seconds. When you see the blinking image with the user number release **⌚** («TIME») button and press it one more time. The blinking year will appear in the display.
- 2) Now You may enter the current year, pressing the **M** («Memory») button. To confirm and set the month press the **⌚** («TIME») button.
- 3) You can now set the month using the **M** («Memory») button. Press the **⌚** («TIME») button to confirm and then set the day.
- 4) Please follow the instructions above to set the day, hour and minutes.
- 5) Once you have set the minutes and pressed the **⌚** («TIME») button, the date and time are set and the time is displayed.

➤ Further Information

With each press of the button (**⌚**), **M**) one input is made (e.g. switching over from hours to minutes mode, or altering the value by +1). However, if you keep the respective button depressed, you can switch more quickly to the desired mode, or find the desired value respectively.

4.5. Specifying user

The WA-55 is designed for two users and can store up to 60 blood pressure readings for each user.

To specify the user, click  («TIME») button and hold it until you see a blinking symbol of the user. To switch between User 1 and User 2 press  («Memory») button. To confirm the choice of the user press Start button.

5. Carrying out a measurement

5.1. Before the measurement:

Avoid eating, smoking as well as all forms of exertion directly before the measurement. All these factors influence the measurement result. Try and find time to relax by sitting in an armchair in a quiet atmosphere for about ten minutes before the measurement.

Remove any garment that fits closely to your upper arm.

Measure always on the same arm (normally left).

Attempt to carry out the measurements regularly at the same time of day, since the blood-pressure changes during the course of the day.

5.2. Common sources of error:

Note:

Comparable blood-pressure measurements always require the same conditions! These are normally always quiet conditions.

- All efforts by the patient to support the arm can increase the blood-pressure. Make sure you are in a comfortable, relaxed position and do not activate any of the muscles in the measurement arm during the measurement. Use a cushion for support if necessary.
- If the arm artery lies considerably lower (higher) than the heart, an erroneously higher (lower) blood-pressure will be measured! (Each 15cm difference in height results in a measurement error of 10mmHg)
- Cuffs that are too narrow or too short result in false measurement values. Selecting the correct cuff is of extraordinary importance. The cuff size is dependent upon the circumference of the arm (measured in the centre).

Note: only use clinically approved Original-Cuffs!

- A loose cuff or a sideways protruding air-pocket causes false measurement values.
- With repeated measurements, blood accumulates in the respective arm, which can lead to false results. Correctly executed blood-pressure measurements should therefore first be repeated after a 1 minute pause.

5.3. Fitting the cuff

- Pass the end of the cuff (with sewn-in rubber stopper) through the metal stirrup so that a loop is formed. The Velcro closer must be facing outwards. (Ignore this step if the cuff has already been prepared.)
- Push the cuff over the left upper arm so that the tube points in the direction of the lower arm.
- Lay the cuff on the arm as illustrated. Make certain that the lower edge of the cuff lies approximately 2 to 3 cm above the elbow and that the rubber tube leaves the cuff on the inner side of the arm. Important! The mark (ca. 3 cm long bar) must lie exactly over the artery which runs down the inner side of the arm.
- Tighten the free end of the cuff and close the cuff with the closer.
- There must be no free space between the arm and the cuff as this would influence the result. Clothing must not restrict the arm. Any piece of clothing which does (e.g. a pullover) must be taken off.
- Secure the cuff with the Velcro closer in such a way that it lies comfortably and is not too tight. Lay the arm on the table (palm upwards) so that the cuff is at the same height as the heart. Make sure that the tube is not kinked.
- Remain seated quietly for two minutes before you begin the measurement.

Comment

If it is not possible to fit the cuff to the left arm, it can also be placed on the right one. However all measurements should be made using the same arm.


5.4. Measurement mode selection

3 Check Mode:

- If you would like to take a 3 Check measurement, please slide the switch to the upside toward no. 3.
- 3Check mode takes generally 3 measurements in succession and calculates the result.
- The « **3CHECK** » symbol in the display indicates that the device is set to the 3 Check mode.

Single mode:

- If you would like to take a single measurement, please slide the switch toward no. 1.
- Single mode only has 1 measurement.

5.5. Measuring procedure

After the cuff has been appropriately positioned the measurement can begin:

- Press the  Start button, pump begins to inflate the cuff. In the display, the increasing cuff-pressure is continually displayed.
- After reaching the inflation pressure, the pump stops and pressure slowly falls away. The cuff-pressure (large characters) is displayed during the measurement.
- When the device has detected the pulse, the heart symbol in the display begins to blink and a peep tone is audible for every pulse beat.
- When the measurement has been concluded, a long peep tone sounds. The measured systolic and diastolic blood-pressure values as well as the pulse frequency are now displayed.
- The appearance of this symbol  signifies that an irregular heartbeat was detected. This indicator is only a caution. It is important that you be relaxed, remain still and do not talk during measurements.
NOTE: We recommend contacting your physician if you see this indicator frequently.
- The measurement results are displayed, until you switch the device off. If no button is pressed for 1 minutes, the device switches automatically off, to save the batteries.
- When the device is set to 3 Check mode, generally 3 separate measurements will take place in succession and calculates your detected blood pressure value. There will be 15 seconds resting time in-between each measurement. A count down indicates the remaining time and a beep will sound 5 seconds before the 2nd and 3rd reading begin. In the case where a single piece of data from each cycle differs too much from each other, a fourth measurement is performed before the result will be displayed. In a rare case in which the blood pressure is unstable even after the fourth measurement <<ERR6>> is shown and no results will be given. If one of the measurements causes an error message, it is repeated.


5.6. Discontinuing a measurement


If it is necessary to interrupt a blood pressure measurement for any reason (e.g. the patient feels unwell), the  Start button can be pressed at any time. The device then immediately lowers the cuff-pressure automatically.

5.7. Memory – storage and recall of the measurements

The blood-pressure monitor automatically stores each of the last 60 measurement values with date & time. By pressing the  «Memory» button, the last measurement (MR60) as well as the further last 59 measurements can be displayed one after the other. If an irregular heartbeat interval was detected in any measurements stored in memory, the icon  will be displayed; Data obtained in 3 Check mode are indicated by a  symbol. To see the measurement results of the second user first switch to the second user (see point 4.5.), then repeat the instructions above.

(M60: Values of the last measurement)

(M59-M1: Values of the measurement before)

5.8. Clear all values

Memory full

When the memory has stored 60 results, the memory is full. From this point onwards, a new measured value is stored by overwriting the oldest value.

Clear all values

If you are sure that you want to permanently remove all stored values, hold down the  «Memory» button (the instrument must have been switched off beforehand) until «CL» appears and then release the button.

To remove all values of the second user first switch  to the second user (see point 4.5.), then repeat the instructions above. Individual values cannot be cleared.

6. Error messages/malfunctions

If an error occurs during a measurement, the measurement is discontinued and a corresponding error code is displayed (Example: Error No. 1).

Error No.	Possible cause(s)
Err 1	No pulse has been detected.
Err 2	Unnatural pressure impulses influence the measurement result. Reason: The arm was moved during the Measurement (Artefact).
Err 3	The inflation of the cuff takes too long. The cuff is not correctly seated.
ERR 5	The measured readings indicated an unacceptable difference between systolic and diastolic pressures. Take another measurement following directions carefully. Contact you doctor if you continue to get unusual readings.
ERR 6	There were too many errors during the measurement in 3Check mode, making it impossible to obtain a final result. Read through the checklist for performing reliable measurements and then repeat the measurement.
«LO»	The pulse is too low (less than 40 beats per minute). Repeat the measurement.
«HI»	The pressure in the cuff is too high (over 300 mmHg) OR the pulse is too high (over 200 beats per minute). Relax for 5 minutes and repeat the measurement.*

Other possible malfunctions and their elimination

If problems occur when using the device, the following points should be checked and if necessary, the corresponding measures are to be taken:

Malfunction	Remedy
The display remains empty when the instrument is switched on although the batteries are in place.	<ol style="list-style-type: none"> 1. Check batteries for correct polarity and if necessary insert correctly 2. If the display is unusual, re-insert batteries or exchange them.
The device frequently fails to measure the blood pressure values, or the values measured are too low (too high).	<ol style="list-style-type: none"> 1. Check the positioning of the cuff. 2. Measure the blood-pressure again in peace and quiet under observance of the details made under point 5.
Every measurement results in different values, although the device functions normally and normal values are displayed.	Please pay attention to the following information as well as the points listed under «Common sources of error» and repeat the measurement.
The blood pressure values measured differ from those measured by the doctor	Record the daily development of the values and consult your doctor about them

► Further Information

The level of blood-pressure is subject to fluctuations even with healthy people. Important thereby is, that comparable **measurements always require the same conditions (Quiet conditions)!**

You must consult your specialist dealer or chemist if there are technical problems with the blood pressure instrument. **Never attempt to repair the instrument yourself! Any unauthorized opening of the instrument invalidates all guarantee claims!**

7. Care and maintenance, recalibration

a) Do not expose the device to either extreme temperatures, humidity, dust or direct sunlight.

b) The cuff contains a sensitive air-tight bubble. Handle this carefully and avoid all types of straining through twisting or buckling.

c) Cleaning the cuff


Clean the device with a soft, dry cloth.

You can machine wash the cuff cover at 30°C (do not iron!).

ⓘ WARNING: Under no circumstances, however, may you wash the inner bladder! Always remove the sensitive bladder from the sleeve before washing and replace it carefully again afterwards.

d) Do not drop the instrument or treat it roughly in any way. Avoid strong vibrations.

e) Never open the device! Otherwise the manufacturer calibration becomes invalid!

Periodical recalibration

Sensitive measuring devices must from time to time be checked for accuracy. We therefore recommend a periodical inspection of the static pressure display every 2 years. Your specialist dealer would be pleased to provide more extensive information about this.

20

8. Reference to standards

Device standard:

Device corresponds to the requirements of the European standard for non-invasive blood- pressure monitor
EN1060-1
EN1060-3
EN1060-4

DIN 58130, NIBP – clinical investigation
ANSI / AAMI SP10, NIBP – requirements

Electromagn. compatibility: Device fulfils the stipulations of the European standard EN 60601-1/-1-2
The stipulations of the EU-Guidelines 93/42/EEC for Medical Products Class IIa have been fulfilled.

9. Technical specifications:

Weight:	580g (with batteries + cuff)
Size:	125(L)x91(W)x59(H)mm
Storage temperature:	-5 to +50°C
Humidity:	15 to 85% relative humidity maximum
Operation temperature:	10 to 40°C
Display:	LCD-Display (Liquid Crystal Display)
Measuring method:	oscillometric
Pressure sensor:	capacitive
Measuring range:	SYS/DIA: 20 to 280 mmHg Pulse: 40 to 200 per minute
Cuff pressure display range:	0-299 mmHg
Memory:	Storing the last 60 measurements for each of 2 users automatically
Measuring resolution:	1 mmHg
Number of independent users:	2

Static Accuracy:	Pressure within ± 3 mmHg
Pulse Accuracy:	± 5 % of the reading
Power source:	4 dry cells (batteries), size AA, 1.5
AC Adapter (optional):	AD-155
Accessories:	EA- Cuff M – L, for arm circumference 22-42 cm or EA-Cuff M, for arm circumference 22-32 cm, washable, storage bag, 4 dry cells (Batteries), size AA, instruction manual, warranty card.

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ONBO Electronic (Shenzhen) Co., Ltd, No. 497, Ta Laneg Nan Road,